




# PUMA 3D Printing Guide

Document created: 26.02.2021, Last edited: 28.06.2025

PUMA: The Open Source (GPL v3.0) 3D printed microscope system designed for direct eye observations and ultra-portability with advanced options for digital imaging, measurements and computer control.



# Contents

Legal Information.....	12
License.....	12
Limitations of Use.....	12
Disclaimer.....	12
Health and Safety Advice.....	13
Burns and Fire.....	13
Toxic Fumes.....	13
Lacerations.....	13
Injury to Eyesight.....	13
Protecting the Vulnerable.....	13
Introduction.....	14
Cover Illustration.....	15
Centre.....	15
From Top Going Clockwise.....	15
Generating Files for 3D Printing from the FreeCAD Source Files.....	16
Generating .3mf files from the models.....	16
Generating STL files from the models.....	17
Printer and Slicer Software and Settings.....	23
Supports.....	23
Mesh Errors.....	24
Post Processing and Aftercare of Prints.....	25
AR Projector.....	26
Resources.....	26
AR_Cc_locknut.....	27
AR_Cc_mount_drawtube.....	28
AR_Cc_Retainer.....	29
AR_Clip.....	30
AR_Light_block_filter.....	31
AR_Cx_Cc_connector.....	32
AR_Cx_Collar.....	33
AR_Cx_Lens_holder.....	34
AR_Stay_thumbwheel.....	35
AR_TFT_DrawTube_F.....	36

AR_TFT_Drawtube_sheath.....	37
AR_TFT_LightShield.....	38
AR_TFT_Mount.....	39
AR_Slide_mount.....	40
Binocular Head.....	41
Resources.....	41
BN_Cap_spacer_bolt_arm.....	42
BN_Cap_spacer_nut_arm.....	43
BN_CM_BB_Cover.....	44
BN_CM_Tube_short_c_adhesin.....	45
BN_EM_Cover.....	46
BN_MB_Shim.....	47
BN_MBT_Cap.....	48
BN_MBT_Ring.....	49
BN_MB_Tube.....	50
BN_MB_Tube_spacer_0p48mm.....	51
BN_MB_Tube_spacer_2p5mm.....	52
BN_Mirror_mount_bottom.....	53
BN_Mirror_mount_top_ocular.....	54
BN_Ocular_spacer_ring_0p32.....	55
BN_Ocular_Base_thread_L.....	56
BN_Ocular_Base_thread_R.....	57
BN_Outlet_lower_spacer.....	58
BN_Outlet_upper.....	59
BN_Splitter_block.....	60
BN_Splitter_support.....	61
Dominus Illumination System.....	62
Location of the models.....	62
Resources.....	65
DI_Cnd_Adj_thumbwheel.....	69
DI_Cnd_Crosshair_Cap_c_adhesin.....	70
DI_Cnd_Flange_Spacer_0p48.....	71
DI_Cnd_gripper.....	72
DI_Cnd_Mirror_holder_socket.....	73
DI_Cnd_protector_cap.....	74

DI_Cnd_to_UC.....	75
DI_Cnd_to_UC_long.....	76
DI_Collector_mirror_block.....	77
DI_Mirror_block_spacer.....	78
DI_Condenser_23_30.....	79
DI_Condenser_base.....	80
DI_Daylight_Kohler_adaptor.....	81
DI_Epi_attachment.....	82
DI_Epi_black_body.....	83
DI_Epi_cap.....	84
DI_Epi_Cnd_Aperture_04.....	85
DI_Epi_Cnd_Aperture_10.....	86
DI_Epi_Cnd_Aperture_13.....	87
DI_Epi_condenser.....	88
DI_Epi_pol.....	89
DI_Ferrule_46_c_adhesin.....	90
DI_Field_Stop_Rectangle.....	91
DI_Field_Stop_Round.....	92
DI_HNA_Diffuser.....	93
DI_HNA_Flange_Spacer_0p48.....	94
DI_HNA_Illuminator.....	95
DI_IAD_AP_25mm_oil.....	96
DI_IAD_AP_DF_Stop_10.....	97
DI_IAD_AP_DF_Stop_11.....	98
DI_IAD_AP_DF_Stop_12.....	99
DI_IAD_AP_DF_Stop_2p5.....	100
DI_IAD_AP_DF_Stop_oil_20.....	101
DI_IAD_AP_Phase180_oil.....	102
DI_IAD_AP_Phase180.....	103
DI_IAD_Filter_Tray_oil.....	104
DI_IAD_Filter_Tray.....	105
DI_IAD_SLM_Cover.....	106
DI_IAD_SLM_Filter.....	107
DI_IAD_Tray.....	108
DI_IFD_Adapter.....	109

DI_IFD_Extension.....	110
DI_IFD_Filter_crosshairs.....	111
DI_IFD_Threadjoiner.....	112
DI_IFD_Threadlock.....	113
DI_IFD_Tray.....	114
DI_LC_Adjust_collar.....	115
DI_LC_Cap.....	116
DI_LC_Collar_2.....	117
DI_LC_Collar_7.....	118
DI_LC_Receptacle.....	119
DI_LC_Spacer.....	120
DI_LED_Cover.....	121
DI_LED_Holder.....	122
DI_LED_Washer.....	123
DI_LPC_Lensholder.....	124
DI_LPC_Lens_retainer.....	125
DI_LPC_Lensless_23.....	126
DI_LPC_Single_23.....	127
DI_M3_Adjustment_ring.....	128
DI_M3_Adjustment_ring_Kohler.....	129
DI_M3_Thumbscrew_short_hex.....	130
DI_M3_Thumbscrew_tall.....	131
DI_Mirror_holder_plain.....	132
DI_Mirror_suspend_plain.....	133
DI_Mirror_to_baseplate.....	134
DI_PWG_CM_Aperture_12.....	135
DI_PWG_Front_stop_12.....	136
DI_PWG_Front_stop_16.....	137
DI_PWG_Pinhole_mount.....	138
DI_PWG_Pinhole_mount_washer.....	139
DI_PWG_Pinhole_mount_well.....	140
DI_PWG_Pinhole_plain.....	141
DI_PWG_Pinhole_presser.....	142
DI_PWG_RMS_to_LC.....	143
DI_PWG_Tripod.....	144

DI_PWG_WindowCap.....	145
DI_Pol_Adjustment_ring.....	146
DI_Pol_middle.....	147
DI_Pol_Receptacle.....	148
DI_Pol_to_baseplate.....	149
DI_Pol_top_bottom.....	150
DI_Proximal_collector_attachment.....	151
DI_Tool_23.....	152
DI_Tool_44.....	153
DI_UC_Retention_Cap.....	154
DI_UC_Spacer.....	155
DI_Uber_pol.....	156
DI_Youngs_frame.....	157
DI_Youngs_slit.....	158
Filterblock.....	159
Resources.....	159
FB_Filter_block_simple.....	160
FB_Filter_collar_compression_tool.....	161
FB_Filter_collar.....	162
FB_Filter_F17_slider.....	163
FB_Filter_F_slider.....	164
FB_Filter_slider.....	165
FB_Filter_Slot_bottom.....	166
FB_Filter_slot_stopper.....	167
FB_Filter_slot_top.....	168
FB_Infinity_adapter.....	169
FB_Side_port_separate.....	170
FB_Side_port_separate_EpiStop.....	171
FB_Splitter_case.....	172
FB_Stopper.....	173
FB_Top_connector.....	174
Focus Gears.....	175
Resources.....	175
FG_Fine_gear.....	176
FG.Focus_spacer.....	177

FG_Intermedius.....	178
FG_Pulley_coarse.....	179
FG_Pulley.....	180
FG_Eccentric_Tensioner_Top.....	181
FG_Eccentric_Tensioner_Bottom.....	182
FG_Eccentric_Tensioner_Pulley_c_adhesin.....	183
Legs.....	184
Resources.....	184
LG_Back_leg_angled.....	185
LG_Feet_linker.....	186
LG_Front_legs.....	187
LG_Hind_extension.....	188
LG_Short_leg.....	189
Monocular Head.....	190
Resources.....	190
MN_EM_Block.....	191
MN_EM_Block_cover.....	192
MN_Monocular_tube_c_adhesin.....	193
MN_Monocular_tube_CM_c_adhesin.....	194
MN_Ocular_cap_170.....	195
MN_Ocular_cap.....	196
MN_Ocular_extension_c_adhesin.....	197
MN_Ocular_extension_CM_c_adhesin.....	198
MN_Ocular_Extn_CM_170_c_adhesin.....	199
MN_Ocular_lock_nut.....	200
MN_Ocular_protective_cap.....	201
MN_Ocular_tube_protective_cap.....	202
MN_Projector_cone.....	203
MN_Aperture_20mm.....	204
MN_Aperture_46mm.....	205
PUMA Control Console.....	206
Resources.....	206
PC_Battery_cover.....	207
PC_Expansion_port_cover.....	208
PC_Front_panel.....	209

PC_Joystick_PCB_clasp.....	210
PC_Lamp_insulator.....	211
PC_Left_panel_Base_Skeleton.....	212
PC_Top_back_panel.....	213
PC_Ardु_cover.....	215
PC_Current_knob.....	216
PC_Monitor_case.....	217
PC_Monitor_connector.....	218
PC_Monitor_lens_retainer.....	219
PC_Monitor_light_shield.....	220
PUMA Lite.....	221
Resources.....	221
PL_Battery_cover.....	222
PL_Knob.....	223
PL_Lamp_insulator.....	223
PL_Left_Base_Back.....	224
PL_Top_Front_Right.....	225
Quick Release Objective Holder.....	227
Resources.....	227
QR_Base_thread.....	228
QR_C-RMS_Thread.....	229
QR_Male_C_extn_1mm.....	230
QR_Male_C_extn_4mm.....	231
QR_Trainer.....	232
QR_Trainer_Male.....	233
Stabiliser.....	234
Resources.....	234
SB_S1_Axial.....	235
SB_S1_Insert_c_adhesin.....	236
SB_S1_Lateral_strut_Left.....	237
SB_S1_Lateral_strut_Right.....	239
SB_S1_Peg.....	240
Stage.....	241
Resources.....	241
ST_Articulation.....	242

ST_BasePlate.....	243
ST_FocusPlate.....	244
Trinocular Camera Port.....	245
Resources.....	245
TN_Trinoc_CP_tube.....	246
TN_Trinoc_CP_CM.....	247
TN_Ocular_spacer_ring_1p44.....	248
XY Stabiliser.....	249
Resources.....	249
XY_Stabiliser.....	249
Z-Motor.....	250
Resources.....	250
ZM_Motor_attachment.....	251
ZM_Motor_gear.....	252
ZM_Z_Limit_Sw_Mount.....	253
ZM_Z_Probe.....	254
ZM_Z_Yoke.....	255
CNC-Stage.....	256
Location of the models.....	256
Resources.....	259
CN_PSU_CableFixture_front.....	263
CN_PSU_DC_CableCover.....	264
CN_PSU_Mains_CableCover.....	265
CN_PSU_Peg_Leg.....	266
CN_Ardu_bus_bar_mount.....	267
CN_Ardu_Bus_Clip_set.....	268
CN_Ardu_M3_Washer.....	269
CN_X_Lim_Sarcoph_L.....	270
CN_X_Lim_Sarcoph_R.....	271
CN_X_NEMA11_Bracket.....	272
CN_X_NEMA11_Counter_weight.....	273
CN_X_NEMA11_CW_Hopper.....	274
CN_X_Probe.....	275
CN_X_Probe_Housing_L.....	276
CN_X_Probe_Housing_R.....	277

CN_X_Probe_shim.....	278
CN_XY_NEMA11_Undertable_Support.....	279
CN_Y_Lim_Sarcoph_Susp_B.....	280
CN_Y_Lim_Sarcoph_Susp_F.....	281
CN_Y_NEMA11_Bracket.....	282
CN_Y_Probe_Susp_Nemesh.....	283
CN_Y_Suspension_platform_B.....	284
CN_Y_Suspension_platform_F.....	284
CN_Z_Coupler.....	285
CN_Z_Lim_Attachment.....	286
CN_Z_Lim_Back_baffle.....	287
CN_Z_Lim_Probe.....	288
CN_Z_M4_small_washer.....	289
CN_Z_Mount.....	290
CN_Z_N11_Syringe_body.....	291
CN_Z_N11_Syringe_cover.....	292
CN_Z_N11_Syringe_Gearbox_plunger.....	293
CN_Z_N11_Syringe_mount.....	294
CN_Sam_Counter_weight.....	295
CN_Sam_Pincer_Extn1_L.....	296
CN_Sam_Pincer_Extn1_R.....	296
CN_Sam_Pincer_Extn2_L.....	297
CN_Sam_Pincer_Extn2_R.....	297
CN_Sam_Presser_wings.....	298
CN_Sam_Spring_arm_assy.....	299
CN_Sam_Stage_Top.....	300
CN_BB_Jig_Back_Left.....	301
CN_BB_Jig_Back_Measure.....	302
CN_BB_Jig_Back_Right.....	303
CN_BB_Jig_Front_Left.....	304
CN_BB_Jig_Front_Measure.....	305
CN_BB_Jig_Front_Right.....	306
CN_BB_Jig_Left_Measure.....	307
CN_BB_Jig_Right_Measure.....	307
CN_BB_Jig_Undercarriage.....	308

CN_BB_Jig_XYTable.....	309
CN_BB_Jig_Y_Lim.....	310
CN_F_ArdU_strut_BL.....	311
CN_F_ArdU_strut_BR.....	311
CN_F_ArdU_strut_TL.....	311
CN_F_ArdU_strut_TR.....	311
CN_F_DM3_Mount.....	312
CN_F_OM_QR_Plate.....	313
CN_F_Overmount.....	314
CN_F_UM_Opt_gripper.....	315
CN_F_UM_Opt_Gripper_washer.....	316
CN_F_UM_Opt_to_QR_F.....	317
CN_F_Undercarriage_platform.....	318
CN_SD_2020_End_cap.....	319
CN_SD_2020_Protector.....	320
CN_SD_Aft_clip.....	321
CN_SD_Aft_Hole_cover.....	321
CN_SD_Aft_Left.....	322
CN_SD_Aft_Mid.....	323
CN_SD_Aft_Right.....	324
CN_SD_Aft_Rivet_block.....	325
CN_SD_DIN_End_Protector.....	326
CN_SD_DIN_Foot_Protector.....	327
CN_SD_Dorsal_Key_head.....	328
CN_SD_Dorsal_Key_key.....	329
CN_SD_Dorsal_Key_main.....	330
CN_SD_Dorsal_Key_wedge.....	330
CN_SD_Dorsal_L.....	331
CN_SD_Dorsal_M.....	332
CN_SD_Dorsal_R.....	333
CN_SD_Dorsal_Rivet_plate.....	334
CN_SD_Dorsal_suture.....	335
CN_SD_Forward_L.....	336
CN_SD_Forward_R.....	336
CN_SD_Forward_Rivet_handle_1.....	337

CN_SD_Forward_Rivet_handle_2.....	337
CN_SD_Forward_Rivet_lock.....	337
CN_SD_Forward_Rivet_tool.....	337
CN_SD_Forward_suture.....	338
CN_SD_Lateral_Left.....	339
CN_SD_Lateral_Right.....	340
CN_SD_Letters_CNC.....	341
CN_SD_Letters_PUMA.....	342
CN_SD_T_nut_positioner.....	343
CN_SD_XMotor_cover_L.....	344
CN_SD_Xmotor_cover_R.....	344

# Legal Information

## License

Copyright (C) 2021-2025 Dr Paul J. Tadrous

Permission is granted to copy, distribute and/or modify this document under the terms of the GNU Free Documentation License, Version 1.3 or any later version published by the Free Software Foundation; with no Invariant Sections, no Front-Cover Texts, and no Back-Cover Texts.

A copy of the license can be found at <https://www.gnu.org/licenses/fdl-1.3.html>

## Limitations of Use

The PUMA microscope and its associated systems do not have any certifications or regulatory approvals in any country for use in clinical diagnostics or treatment (human or veterinary).

The PUMA microscope and its associated systems are released to be used for research and educational purposes only.

## Disclaimer

All PUMA project information, including without limitation any CAD file (or 3-dimensional (3D) model or mesh file in any other format) and all documentation, advice and instruction (whether provided in video form, audible form, written form or otherwise) is provided 'as is' in good faith and is intended to be helpful but comes with no warranty whatsoever.

Anyone attempting to build or use a PUMA microscope or other PUMA-related material, accessory, module or derivative is hereby advised that there will be risk involved in 3D printing, post-print processing, assembly and usage of the resulting structures. This risk includes, without limitation, the risk of personal damage and loss of resources.

Dr Paul J. Tadrous, TadPath and OptArc cannot accept any liability for any such loss or damages that may occur. All those who attempt to build or use any aspect of the PUMA project or derivatives thereof do so at their own risk.

# **Health and Safety Advice**

## **Burns and Fire**

There is a risk of fire and of burns when working with 3D printers, soldering irons and molten solder. Wear appropriate personal protective equipment such as gloves and eye and face protection. Ensure all electrical components are of the appropriate voltage and current tolerances. Do not leave electrical components powered for longer than they are required for use. Do not leave 3D printers unattended when they are powered up.

## **Toxic Fumes**

Some filaments used in 3D printing can release toxic fumes when heated (such as during the 3D printing process). Working with solder can also release toxic fumes (especially if leaded solder is used). Ensure adequate ventilation at all times during 3D printing and soldering.

## **Lacerations**

Cleaning and assembling 3D prints and other components of PUMA systems sometimes involves the use of sharp instruments – take appropriate precautions when handling sharp blades. Working with optics – especially thin plates of glass such as beam splitter mirrors, microscope slides and coverslips and display screen panels can result in flying shards and sharp edges. Wear appropriate personal protective equipment such as gloves and face and eye protection.

## **Injury to Eyesight**

Physical injuries to the eyes may occur from shards and specs of glass and airborne particles of molten or hot solder. Bright lights are used to illuminate the PUMA microscope. These should be used with caution. Do not stare directly into the illumination LED bulb (with or without its collector lenses in place), and only use a brightness setting that is appropriate for comfortable viewing when using the complete microscope. If using a mirror to provide illumination always use a reflection from a diffuse scene such as an open blue or cloudy sky. NEVER allow specular reflections of the sun to enter the field of view of the mirror or permanent eyesight loss may result. If using a bulb that emits UV light, always ensure an appropriate and effective UV blocking filter is used to shield your eyes from the UV light. If a laser is used for illumination then you should never look at the image through the lenses of the scope (you should use a camera for making observations).

## **Protecting the Vulnerable**

The PUMA microscope and its accessories are not toys. Do not allow children, vulnerable adults or pets to come into contact with any PUMA product without appropriate supervision and safeguards.

# Introduction

PUMA is a 3D printed high quality microscope system. The name PUMA stands for **P**ortable, **U**pgradeable, **M**odular and **A**ffordable – key features of this system.

This document provides practical advice for those who want to print their own parts for the building of PUMA-related modules and microscopes.

It also serves as a visual ‘index of parts’ for all the 3D printed parts in the PUMA project and shows which FreeCAD file contains which parts. Within the FreeCAD file those models have the same name as the printable model file but without the printable file extension (like ‘.stl’ or ‘.3mf’) and also without the two letter prefix before the first underscore. For example, the printable file ‘AR\_Cc\_locknut.stl’ is a mesh generated from the model called ‘Cc\_locknut’ in the FreeCAD file called ‘AR\_Projector.FCStd’.

The information in this file only relates to 3D printed parts. To build a PUMA microscope additional components will be needed such as nuts and bolts and optics elements. None of these additional components are addressed here. See the assembly / usage tutorials, videos and other ‘how to build’ documentation for that information.

At the beginning of each section some general printing advice (like slicer settings profile) is given if the models for that section are to be printed with a profile other than the default (to be discussed later). Further advice is given for each model separately. For most models this advice consists solely in an illustration showing the preferred orientation of the model on the printer build plate. In some cases additional advice and illustrations are given, particularly where they relate to support structures.

This document makes no attempt at explaining what each part does or how to use or assemble them. For that kind of information you would need to see the assembly / usage tutorials and other documentation.

The Cura screenshot illustrations in this document show the models at various zooms and with some angular perspective variations from model to model in order to emphasise different aspects. However the grid and ‘Ender’ logo shown on the build plate should make it clear how the model is to be orientated and indicate the scale.

Each chapter begins with a table of resources showing the print time and length of 1.75 mm diameter PLA filament used to print each model file as estimated by Cura v.4.8.0 using the recommended print settings described in that chapter for each model.

# Cover Illustration

## Centre

A PUMA microscope in the following configuration: Full Kohler illuminator, motorised Z-stage, mechanical XY Vernier caliper slide holder, advanced filter block, augmented reality projector, binocular head, PUMA control console. It is fitted with a Zeiss Plan 2.5x objective and two high eye point wide field 10x oculars of 20 mm field of view. Many alternative configurations are possible due to PUMA's modular design.

## From Top Going Clockwise

These microscope images (all taken with a PUMA microscope in various configurations) are highly cropped and compressed and do not show the full field of view. See the gallery of images on the PUMA websites for full quality and full field images. Magnifications given below are of the objective only:

Bright field Kohler illumination: H&E chronic inflammation, x100 oil immersion.

Fluorescence microscopy: Fluorescein stained *Soncus* bud showing pollen grains and other cellular structures x40.

Fluorescence microscopy: Fluorescein stained //ex (holly) leaf showing chloroplasts, vascular structures and other cellular structures x40.


Crossed polarisation microscopy: Amyloid stained with Congo red x40.

Dark field / Dark ground microscopy: Unstained section of colonic tissue x10.

Schlieren phase contrast using a spatial light modulator in the illuminated aperture plane of the condenser: Unstained section of colonic mucosa, x100 oil immersion.

Crossed polarisation microscopy: Striated muscle of tongue stained with H&E x40.

AR HUD Overlay. The figure shows a cropped area from a small part of the field of a x10 view of an H&E section (of testis) with the augmented reality (AR) projector projecting a heads up display (HUD) of a grid and some microscope status information onto the live optical image. The full field image is shown here:




# Generating Files for 3D Printing from the FreeCAD Source Files

In previous editions of this document I used mainly STL mesh files for 3D printing. I have since moved away from favouring STL in favour of 3mf files as I found that STL files could sometimes be excessively large (in terms of megabytes) and might not load properly (or take excessively long to load) in the slicer program that I use (Cura) as well as other slicers. Furthermore, STL file generation with FreeCAD requires an intermediate processing step on the models (to turn them into meshes) but FreeCAD models can be saved directly as 3mf files without any processing. This latter distinction became especially important since some bugs have crept into later models of FreeCAD which do not allow the reading of older FreeCAD files properly once they are imported into the new FreeCAD version and 'recalculated' for processing. This means that, until such bugs are fixed, some of the older models cannot be edited in FreeCAD (such as to make them into meshes) but they can be read and directly saved as 3mf files. For all those reasons I tend to use 3mf files for printing rather than stl files. However, I will show you how to save models in both formats so you can choose (the latest FreeCAD files in the PUMA project – such as those for the CNC stage – are made in the modern (post v.1.0) FreeCAD so may be edited and saved as STL if desired).

## Generating .3mf files from the models

Due to issues with FreeCAD v.1.0 onwards, take the following steps to avoid printing faulty models if using many of the older PUMA models. These steps will also work with the newer PUMA FreeCAD models:


Load the FreeCAD file but do NOT 'recompute' the model as suggested by FreeCAD (and make sure your FreeCAD does not automatically re-compute models on loading). Instead just select the model and export it as a .3mf file by selecting File -> Export (see the figure).




## Generating STL files from the models

If you prefer or need to generate STL files instead of 3mf files, you can do so by means of the following steps. However, these will require you to have fully computed models so you may find this option is not available for the older PUMA models (prior to 2025) due to modern FreeCAD incompatibilities when loading older FreeCAD files. In some cases the files may undergo the 'recompute' process without crashing but, in that case, some aspects of the models may be altered making them 'wrong'. This is why the 3mf file route discussed above is preferred.


1. Load the CAD file into FreeCAD (I used v. 0.18 to illustrate these steps but a later or the latest version may be required to correctly view and edit the latest models). For the following illustrated examples I have loaded the 'Z\_Motor.FCStd' file from the PUMA project.
2. Select the 'Mesh Design' Workbench from the drop down selection list on the top button bar of FreeCAD.




3. Left-click on the model you want to create a mesh for so as to select it.



4. Go to 'Meshes' -> 'Create mesh from shape ...' to bring up the Tessellation dialogue box.




5. Select the 'Standard' meshing option set the 'surface deviation' to 0.001 mm and the 'Angular deviation' to 1 degree.



6. Click OK - it may take a while to complete.

These mesh settings may be overly fine for many of the models in that you could probably get away with coarser settings and lower resolution meshes for some models. However I have not experimented to see how low I can go before the models cease to function or cease to fit properly.

7. When the mesh is generated, right click on the mesh in the 'Model' tab of the combo view and select 'Export Mesh ...' from the drop down menu to save the mesh as an STL file. You might need to rotate the mesh first before export to make it easier to site on the build plate in Cura in the manner required for printing (as will be shown in the illustrations in the rest of this document for each model). Although you can also perform most of those manipulations in Cura itself.



# Printer and Slicer Software and Settings

This advice relates specifically to printing parts using a Creality Ender 3 printer with ordinary PLA 1.75 mm diameter filament and the Ultimaker Cura slicer software (v. 4.8.0 to v. 4.13.1).

If you use a different printer and / or slicer you would need to explore for yourself what are the most appropriate settings and profiles to use and may also need to modify other aspects of the advice given here.

Note the orientation of models on the print bed as shown in the figures. For some models is probably won't matter if you use a different orientation but for others is it crucial to functionality that you stick to the orientation shown. My advice therefore is to always use the orientation shown in the figures unless you want to experiment and are happy to risk either failure or suboptimal functionality.

All models are printed with the PUMA custom Cura profile called 'Flats' unless otherwise specified (several other PUMA custom Cura profiles are used on occasion and these are available in the GitHub repository). No supports are used unless specified otherwise. None of the Cura in-built build plate adhesion options are used. Where necessary, custom adhesion structures – called 'adhesins' – are built into the model and these will need to be trimmed off when the model is complete as part of model post print clean up.

I have only ever printed the parts of this project using super high quality print settings and the highest resolution tessellation meshes. No doubt this is 'overkill' for many of the parts but I have not done the practical experimentation required to find out the lowest resolution meshes that will work with the lowest quality / fastest print settings.

I have only ever used standard PLA 1.75 mm filament for printing the parts of this project – although I have used PLA from different suppliers from Creality's own brand to much cheaper generic filament. With the highest quality branded filaments the prints seem a bit smoother but from a functionality point of view all filaments I have tried give acceptable results.

## Supports

You don't need supports for screw thread holes up to M4 diameter that are printed horizontally but Cura will try to insert supports for these which is one of the reasons why support advice in some of the models below specifies to use the 'touching build plate only' option (not 'Everywhere').

When supports are enabled the default overhang angle should be set to 67 degrees.

For 'Tree' supports the default parameters to use are: 'branch distance' is 1 mm, 'branch diameter' is 2 mm, 'branch diameter angle' is 5 degrees and 'collision resolution' is 0.2 mm.

The 'branch angle' for Tree supports may be 40 or 50 degrees depending on the model so this information is provided for each model that required tree supports.

If supports are not stated to be required then ensure supports are off – even if it appears (or Cura suggests) that they might be needed. Sometimes this is emphasised

for models where having supports enabled could cause the print to fail or adversely affect the function of the resulting print.

Some models require specific support blockers. These will be illustrated in the figures.

## Mesh Errors

Very occasionally for some models Cura will report model errors in the STL mesh and highlight those areas in bright colours. These ‘errors’ are very small and the model slices and prints just fine even with them so they can be ignored. During development I have been able to identify and fix some of those mesh errors in FreeCAD but the very few remaining ones were not obvious – and practically they do not matter for printing with an Ender 3 printer.

## Post Processing and Aftercare of Prints

I have noticed that PLA prints will shrink a little if exposed to excessive heat post printing and this may distort closely fitting parts and prevent them from functioning adequately. To avoid this, never keep your 3D printed parts on or very near to a radiator and never use a heater, hair dryer or radiator to dry off any parts that may have gotten wet. Leaving parts in a car in the sun may have a similar detrimental heat shrinkage effect.

Some parts, especially closely fitting parts like threads and push fit components may need scraping or light sanding to remove 'nerds', 'zits' or prominent seams or other imperfections prior to use. It is advised that you allow parts to cool down thoroughly after printing before any such actions (especially sanding) otherwise you may find that you have over-eroded the model and end up with a fit that is too loose. Some 3D printed threads may require 'training' if they seem overly tight. This means you screw and unscrew the thread repeatedly to wear it down a bit – this will usually wear down small prominences / defects in the print that lie in the thread channels. This may need to be done gradually (you may not be able to screw a part all the way in at first). After training the thread should fit snugly without being too tight or too loose.

Some parts require lubrication – such as the polariser wheel and Z-stage bearings. This will be elaborated on in the construction and usage tutorials.

Many thread holes for screws (mostly M2 and M3 sizes but sometimes also up to M6) in many of the parts are printed simply as holes that are slightly too small for the screw and you are expected to cut a thread into the hole by simply screwing the metal screw into the hole. Do this slowly at first concentration on getting the screw perpendicular to the hole and applying downward pressure. For such joints it is important never to overtighten the screws in the finished product because you can easily destroy the plastic thread so produced. Also, you should avoid unscrewing and re-screwing in those holes because the thread will wear out eventually (although they can take quite a few repeats before wearing down to a detrimental degree).

## AR Projector


These are the parts for the augmented reality (AR) projector module which provides a heads-up-display (HUD) graphical user interface to the user. The CAD source models for these files are found in the file AR\_Projector.FCStd.

## Resources


Cura calculates the following resources are required to print each model in this chapter:

AR_Projector	Time_Hr	Time_Min	PLA_Length(m)
AR_Cc_locknut	0	12	0.35
AR_Cc_mount_drawtube	1	49	6.3
AR_Cc_Retainer	0	7	0.28
AR_Cx_Cc_connector	0	28	1.16
AR_Cx_Collar	0	11	0.36
AR_Cx_Lens_holder	0	32	1.97
AR_Clip	0	45	1.41
AR_Light_block_filter	0	1	0.06
AR_Stay_thumbwheel	0	15	0.56
AR_TFT_DrawTube_F	2	59	5.33
AR_TFT_Drawtube_sheath	0	18	0.87
AR_TFT_LightShield	0	43	1.49
AR_TFT_Mount	1	1	2.3
AR_Slide_mount	0	59	1.91

## AR\_Cc\_locknut




## AR\_Cc\_mount\_drawtube




Print with PUMA custom Cura profile 'P\_AR\_Tube'

## AR\_Cc\_Retainer




## AR\_Clip




Use standard supports (not tree) with overhang of 67 degrees.


## AR\_Light\_block\_filter




## AR\_Cx\_Cc\_connector



## AR\_Cx\_Collar




## AR\_Cx\_Lens\_holder




Print with PUMA custom Cura profile 'Strong01'

## AR\_Stay\_thumbwheel




Use profile 'Strong01'. No supports.


## AR\_TFT\_Drawtube\_F




## AR\_TFT\_Drawtube\_sheath




## AR\_TFT\_LightShield



## AR\_TFT\_Mount



## AR\_Slide\_mount



## Binocular Head


These are the 3D printed parts for the binocular head module. The CAD source models for these files are found in the file Binocular.FCStd.

## Resources


Cura calculates the following resources are required to print each model in this chapter:

Binocular_Head	Time_Hr	Time_Min	PLA_Length(m)
BN_Cap_spacer_bolt_arm	0	21	0.77
BN_Cap_spacer_nut_arm	0	22	0.95
BN_CM_BB_Cover	0	31	1.13
BN_CM_Tube_short_c_adhesin	0	41	1.69
BN_EM_Cover	0	21	0.87
BN_MB_Shim	0	5	0.08
BN_MBT_Cap	0	24	1
BN_MBT_Ring	0	21	0.64
BN_MB_Tube	1	50	3.68
BN_MB_Tube_spacer_2p5mm	0	12	0.28
BN_Mirror_mount_bottom	1	0	1.51
BN_Mirror_mount_top_ocular	1	31	2.84
BN_Ocular_spacer_ring_0p32	0	2	0.06
BN_Ocular_Base_thread_L	0	35	0.86
BN_Ocular_Base_thread_R	0	35	0.86
BN_Outlet_lower_spacer	0	0	0.02
BN_Outlet_upper	1	42	2.72
BN_Splitter_block	2	4	3.35
BN_Splitter_support	0	1	0.02

## BN\_Cap\_spacer\_bolt\_arm




## BN\_Cap\_spacer\_nut\_arm




Use the Condenser profile (including its support settings).


## BN\_CM\_BB\_Cover




## BN\_CM\_Tube\_short\_c\_adhesin




## BN\_EM\_Cover




## BN\_MB\_Shim




## BN\_MBT\_Cap




## BN\_MBT\_Ring




## BN\_MB\_Tube



## BN\_MB\_Tube\_spacer\_0p48mm




## BN\_MB\_Tube\_spacer\_2p5mm




This is just one example from a range of MB\_Tube\_spacers of various thicknesses.

## BN\_Mirror\_mount\_bottom




Select supports 'touching build plate only' (not 'Everywhere').

## BN\_Mirror\_mount\_top\_ocular




Select supports 'touching build plate only' (not 'Everywhere').

## BN\_Ocular\_spacer\_ring\_0p32




Use 'Top/bottom pattern'='Concentric' in the Top/Bottom' settings

## BN\_Ocular\_Base\_thread\_L




Select supports 'touching build plate only' (not 'Everywhere').

## BN\_Ocular\_Base\_thread\_R




Select supports 'touching build plate only' (not 'Everywhere').

## BN\_Outlet\_lower\_spacer




## BN\_Outlet\_upper




Use supports of 'Tree' type 'touching build plate only' (not 'Everywhere') and change the default for 'tree support branch distance' from 1.0 mm to 0.2 mm.

## BN\_Splitter\_block



Select supports 'touching build plate only' (not 'Everywhere').

## BN\_Splitter\_support



# Dominus Illumination System

These are parts for the modular illumination system of the microscope.

The illumination system includes an extensive range of modules for various types of illumination and has many more models than other aspects of the PUMA microscopy system. If all placed in 1 file it would go over the 25 MB limit for files on GitHub so these models are split into 3 separate FreeCAD files: Dominus\_part1.FCStd, Dominus\_part2.FCStd and Dominus\_part3.FCStd.

Also, because there are so many models, to make it easier for users to identify which file a particular model is in I have put the models into separate functionally related groups within each file and a list of where each model can be found is given below. Note that some modules require models from different functional groups and the instructional videos and other 'How To' documentation will make it clear what models are needed to build any particular illumination module.

## Location of the models

### File: Dominus\_part1.FCStd

#### *Group: LED\_Housing*

- LED\_Washer
- LED\_Cover
- LED\_Holder

#### *Group: Mirror\_illumination*

- Mirror\_suspend\_plain
- Mirror\_to\_baseplate
- Mirror\_holder\_plain
- Cnd\_mirror\_holder\_socket
- Daylight\_Kohler\_adaptor

#### *Group: IFD*

- IFD\_Threadlock
- IFD\_Tray
- Field\_Stop\_Round
- Field\_Stop\_Rectangle
- IFD\_Filter\_Crosshairs
- IFD\_Extension\_c\_adhesin
- M3\_Adjustment\_ring\_Kohler

*Group: Lower\_Collector*

- LC\_Cap
- LC\_Receptacle
- LC\_Adjust\_collar
- LC\_Collar\_2
- LC\_Collar\_7
- LC\_Spacer

*Group: Tools*

- Tool\_23
- Tool\_44

*Group: Low\_Power\_Collimator*

- IFD\_Adapter
- IFD\_Threadjoiner
- Ferrule\_46\_c\_adhesin
- LPC\_Single\_23
- LPC\_Lensholder
- LPC\_Lens\_retainer
- LPC\_Lenseless\_23

*Group: M3\_Thumbscrews*

- M3\_Thumbscrew\_short\_hex

*Sub-group: M3\_Thumbscrew\_tall*

- M3\_knob

**File: Dominus\_part2.FCStd**

*Group: IAD*

- IAD\_SLM\_Cover
- IAD\_SLM\_Filter
- IAD\_AP\_DF\_Stop\_oil\_20
- IAD\_AP\_Phase180\_oil
- IAD\_Filter\_tray\_oil
- IAD\_AP\_DF\_Stop\_2p5
- IAD\_AP\_DF\_Stop\_11
- IAD\_AP\_DF\_Stop\_10
- IAD\_AP\_DF\_Stop\_12
- IAD\_AP\_Phase180

IAD\_Filter\_Tray

IAD\_Tray

IAD\_Filter\_template

IAD\_AP\_25mm\_oil

*Group: Upper\_Collector*

UC\_Spacer

UC\_Retention\_Cap

Cnd\_to\_UC\_long

Cnd\_to\_UC

*Group: Intercollector\_Mirror\_block*

Mirror\_block\_spacer

Proximal\_collector\_attachment

Collector\_mirror\_block

*Group: Polariser\_Condenser*

Pol\_to\_baseplate

Pol\_Receptacle

Pol\_top\_bottom

Pol\_middle

Pol\_Adjustment\_ring

Uber\_pol

*Group: Epi\_Illumination*

Epi\_Cnd\_Aperture\_10

Epi\_Cnd\_Aperture\_04

Epi\_Cnd\_Aperture\_13

Epi\_condenser

Epi\_condenser

Epi\_attachment

M3\_Adjustment\_ring

Epi\_black\_body

Epi\_pol

*Group: Abbe\_Condenser*

Cnd\_Adj\_thumbwheel

Cnd\_protector\_cap

Cnd\_Crosshair\_Cap\_c\_adhesin

Cnd\_gripper

Cnd\_Flange\_Spacer\_0p48  
Condenser\_base  
Condenser\_23\_30

#### File: Dominus\_part3.FCStd

*Group: HNA\_Illuminator001*

HNA\_Illuminator  
HNA\_Diffuser  
HNA\_Flange\_Spacer\_0p48

*Group: Plane\_Wave\_Generator*

PWG\_Pinhole\_plain  
PWG\_Pinhole\_mount  
PWG\_Tripod  
PWG\_CM\_Aperture\_12  
PWG\_Pinhole\_mount\_washer  
PWG\_Pinhole\_presser\_c\_adhesin  
PWG\_Front\_stop\_12  
PWG\_Front\_stop\_16  
PWG\_WindowCap\_c\_adhesin  
PWG\_RMS\_to\_LC\_c\_adhesin  
PWG\_Pinhole\_mount\_well\_c\_adhesin

*Group: Youngs\_Slits*

Youngs\_slit  
Youngs\_frame

## Resources

Cura calculates the following resources are required to print each model in this chapter:


Dominus_Illuminator	Time_Hr	Time_Min	PLA_Length(m)
DI_Cnd_Adj_thumbwheel	0	8	0.26
DI_Cnd_Crosshair_Cap_c_adhesin	0	21	0.87
DI_Cnd_Flange_Spacer_0p48	0	5	0.19
DI_Cnd_gripper	1	30	3.69
DI_Cnd_Mirror_holder_socket	1	33	4.43

DI_Cnd_protector_cap	0	22	0.98
DI_Cnd_to_UC	0	47	2.89
DI_Cnd_to_UC_long	1	8	3.87
DI_Collector_mirror_block	5	9	12.51
DI_Condenser_23_30	2	31	5.94
DI_Condenser_base	0	46	1.38
DI_Daylight_Kohler_adaptor	1	17	3.25
DI_Epi_attachment	0	47	1.84
DI_Epi_black_body	1	33	4.14
DI_Epi_cap	0	13	0.58
DI_Epi_Cnd_Aperture_04	0	4	0.15
DI_Epi_Cnd_Aperture_10	0	3	0.1
DI_Epi_Cnd_Aperture_13	0	2	0.07
DI_Epi_condenser	2	30	6.23
DI_Epi_pol	0	58	2.02
DI_Ferrule_46_c_adhesin	1	16	4.08
DI_Field_Stop_Rectangle	0	3	0.12
DI_Field_Stop_Round	0	3	0.11
DI_HNA_Diffuser	0	24	0.93
DI_HNA_Flange_Spacer_0p48	0	4	0.14
DI_HNA_Illuminator	0	58	3.08
DI_IAD_AP_25mm_oil	0	1	0.05
DI_IAD_AP_DF_Stop_10	0	1	0.04
DI_IAD_AP_DF_Stop_11	0	1	0.04
DI_IAD_AP_DF_Stop_12	0	1	0.05
DI_IAD_AP_DF_Stop_2p5	0	1	0.03
DI_IAD_AP_DF_Stop_oil_20	0	2	0.08
DI_IAD_AP_Phase180_oil	0	2	0.08
DI_IAD_AP_Phase180	0	1	0.06
DI_IAD_Filter_Tray_oil	0	33	1.41
DI_IAD_Filter_Tray	0	37	1.64
DI_IAD_SLM_Cover	0	19	0.76
DI_IAD_SLM_Filter	0	29	0.93


DI_IAD_Tray	0	27	0.72
DI_IFD_Adapter	1	5	2.88
DI_IFD_Extension_c_adhesin	2	10	7.28
DI_IFD_Filter_crosshairs	0	37	1.59
DI_IFD_Threadjoiner	0	50	1.8
DI_IFD_Threadlock	0	24	0.84
DI_IFD_Tray	2	2	4
DI_LC_Adjust_collar	0	20	0.88
DI_LC_Cap	0	21	1.02
DI_LC_Collar_2	0	2	0.07
DI_LC_Collar_7	0	7	0.22
DI_LC_Receptacle	0	21	0.77
DI_LC_Spacer	0	16	0.44
DI_LED_Cover	0	34	0.87
DI_LED_Holder	0	57	1.87
DI_LED_Washer	0	2	0.06
DI_LPC_Lensholder	0	16	0.54
DI_LPC_Lens_retainer	0	10	0.38
DI_LPC_Lenseless_23	0	31	1.08
DI_LPC_Single_23	0	22	0.81
DI_M3_Adjustment_ring	1	7	2.35
DI_M3_Adjustment_ring_Kohler	2	11	5.2
DI_M3_Thumbscrew_short_hex	0	8	0.14
DI_M3_Thumbscrew_tall	0	19	0.29
DI_Mirror_block_spacer	0	52	1.73
DI_Mirror_holder_plain	1	0	2.45
DI_Mirror_suspend_plain	1	34	3.18
DI_Mirror_to_baseplate	0	49	1.88
DI_Pol_Adjustment_ring	1	22	2.58
DI_Pol_middle	0	5	0.16
DI_Pol_Receptacle	0	46	1.28
DI_Pol_to_baseplate	2	19	4.88
DI_Pol_top_bottom	0	7	0.25

DI_Proximal_collector_attachment	0	38	1.32
DI_PWG_CM_Aperture_12	0	2	0.06
DI_PWG_Pinhole_plain	0	5	0.24
DI_PWG_Pinhole_mount	0	5	0.24
DI_PWG_Pinhole_mount_washer	0	5	0.15
DI_PWG_Front_stop_16	0	18	0.91
DI_PWG_Front_stop_12	0	20	1
DI_PWG_Pinhole_mount_well	0	41	1.43
DI_PWG_Pinhole_presser	0	25	1.17
DI_PWG_RMS_to_LC	0	29	1.35
DI_PWG_Tripod	0	58	2.14
DI_PWG_WindowCap	0	14	0.6
DI_Tool_23	0	24	1.14
DI_Tool_44	0	28	0.84
DI_UC_Retention_Cap	0	10	0.59
DI_UC_Spacer	0	1	0.05
DI_Uber_pol	0	12	0.19


## DI\_Cnd\_Adj\_thumbwheel




## DI\_Cnd\_Crosshair\_Cap\_c\_adhesin




## DI\_Cnd\_Flange\_Spacer\_0p48



## DI\_Cnd\_gripper




## DI\_Cnd\_Mirror\_holder\_socket




This uses tree support touching baseplate only with tree support branch angle 40 degrees and tree support branch distance 1 mm.


## DI\_Cnd\_protector\_cap




## DI\_Cnd\_to\_UC




## DI\_Cnd\_to\_UC\_long




## DI\_Collector\_mirror\_block




Supports on for 'Everywhere' and ensure 'Support overhang angle' is 67. Make support blockers for the M3 holes in the front flat part, the M3 adjustment ring holes and the part of the mirror indentation that overhangs - essentially everywhere except the big vertical aperture - that is what needs support.



## DI\_Mirror\_block\_spacer




## DI\_Condenser\_23\_30




Must be printed with the special custom Cura profile called 'Condenser' which I made to ensure both strength and proper overhang tree supports to the centration flange disc.

## DI\_Condenser\_base




Ensure the 'Top/bottom Pattern' is set to 'Concentric'. This avoids nozzle blebs on the diaphragm surface (they can cause deformations upwards when the condenser attachment is screwed in and such deformations can hinder passage of IAD filters).

## DI\_Daylight\_Kohler\_adaptor




Ensure the 'Top/bottom Pattern' is set to 'Concentric'. I suggest also that the infill pattern be set to 'zig-zag' to reduce retractions.

## DI\_Epi\_attachment




## DI\_Epi\_black\_body




Use the standard 'Flats' profile and do not use supports.


## DI\_Epi\_cap




## DI\_Epi\_Cnd\_Aperture\_04




## DI\_Epi\_Cnd\_Aperture\_10



## DI\_Epi\_Cnd\_Aperture\_13




## DI\_Epi\_condenser




Use the 'Pol\_to\_baseplate' Cura profile.


## DI\_Epi\_pol



## DI\_Ferrule\_46\_c\_adhesin



## DI\_Field\_Stop\_Rectangle




Note that the Field\_Stop models are thicker (have more print layers) than the DI\_IFD\_Filter\_Crosshairs. This is because it is beneficial for the crosshairs to be semitransparent to be able to see the image of the light source behind it for centration purposes but true field stops should be fully opaque.

In FreeCAD:

To change the aperture of the Field\_Stop\_Rectangle alter the dimensions (and XY offset positions accordingly to maintain central position) of the box called 'Aperture\_rectangle' here (the default size is 4 mm by 3 mm):

Field\_Stop\_Rectangle  
  |\_ Sehfeldblende\_4x3  
    |\_ Aperture\_rectangle

## DI\_Field\_Stop\_Round




Note that the Field\_Stop models are thicker (have more print layers) than the DI\_IFD\_Filter\_Crosshairs. This is because it is beneficial for the crosshairs to be semitransparent to be able to see the image of the light source behind it for centration purposes but true field stops should be fully opaque.

In FreeCAD:


To change the aperture of the Field\_Stop\_Round alter the radius value of the cylinder called 'Aperture\_hole' here (the default radius is 3.5 mm to give a 7 mm diameter aperture):

```
Field_Stop_Round
|_ Sehfeldblende
    |_ Aperture_hole
```


## DI\_HNA\_Diffuser



## DI\_HNA\_Flange\_Spacer\_0p48




## DI\_HNA\_Illuminator




Supports on for 'touching build plate only'


## DI\_IAD\_AP\_25mm\_oil




## DI\_IAD\_AP\_DF\_Stop\_10




## DI\_IAD\_AP\_DF\_Stop\_11




## DI\_IAD\_AP\_DF\_Stop\_12




## DI\_IAD\_AP\_DF\_Stop\_2p5




## DI\_IAD\_AP\_DF\_Stop\_oil\_20



## DI\_IAD\_AP\_Phase180\_oil



## DI\_IAD\_AP\_Phase180




## DI\_IAD\_Filter\_Tray\_oil




## DI\_IAD\_Filter\_Tray




## DI\_IAD\_SLM\_Cover




## DI\_IAD\_SLM\_Filter




## DI\_IAD\_Tray




## DI\_IFD\_Adapter




## DI\_IFD\_Extension




## DI\_IFD\_Filter\_crosshairs




## DI\_IFD\_Threadjoiner



## DI\_IFD\_Threadlock




## DI\_IFD\_Tray




Supports on for 'touching build plate only'


## DI\_LC\_Adjust\_collar




## DI\_LC\_Cap




## DI\_LC\_Collar\_2




## DI\_LC\_Collar\_7




## DI\_LC\_Receptacle




## DI\_LC\_Spacer




## DI\_LED\_Cover




## DI\_LED\_Holder




## DI\_LED\_Washer




## DI\_LPC\_Lensholder




## DI\_LPC\_Lens\_retainer




## DI\_LPC\_Lensless\_23



## DI\_LPC\_Single\_23




## DI\_M3\_Adjustment\_ring



Ensure all supports are OFF


## DI\_M3\_Adjustment\_ring\_Kohler




Ensure all supports are OFF.

Ignore the 'model error' / 'mesh error' warnings.


## DI\_M3\_Thumbscrew\_short\_hex




## DI\_M3\_Thumbscrew\_tall



## DI\_Mirror\_holder\_plain




## DI\_Mirror\_suspend\_plain




Use the 'Plates' custom Cura profile and enable supports 'touching baseplate only'.

## DI\_Mirror\_to\_baseplate




## DI\_PWG\_CM\_Aperture\_12




Use 'WasherSpacer' Cura profile.


## DI\_PWG\_Front\_stop\_12




## DI\_PWG\_Front\_stop\_16




## DI\_PWG\_Pinhole\_mount




## DI\_PWG\_Pinhole\_mount\_washer




## DI\_PWG\_Pinhole\_mount\_well




## DI\_PWG\_Pinhole\_plain




## DI\_PWG\_Pinhole\_presser




## DI\_PWG\_RMS\_to\_LC




## DI\_PWG\_Tripod




## DI\_PWG\_WindowCap




## DI\_Pol\_Adjustment\_ring




## DI\_Pol\_middle



## DI\_Pol\_Receptacle




## DI\_Pol\_to\_baseplate




Pol\_to\_baseplate has its own Cura profile with tree supports. This is similar to the 'Condenser' profile but uses an infill density of 20% (not 30%) and has a greater tree angle to avoid making supports on the base of the model.


## DI\_Pol\_top\_bottom



## DI\_Proximal\_collector\_attachment




## DI\_Tool\_23




Although shown here with profile 'flats' you may want to print this with a higher density infill for added strength.


## DI\_Tool\_44



## DI\_UC\_Retention\_Cap




## DI\_UC\_Spacer




You may find the custom Cura profile 'WashersSpacers' more useful when printing this.


## DI\_Uber\_pol



## DI\_Youngs\_frame



## DI\_Youngs\_slit



Note the orientation of the part is important, not only to ensure a good straight optical edge by motion along the X-gantry only but also to avoid Z-seam lines at the optical flat edge.

# Filterblock


This shows the parts for the filter block segment of the optical tube, the part that goes between the objective holder and the ocular head. The models for these files are found in the FreeCAD file called FilterBlock.FCStd.

## Resources


Cura calculates the following resources are required to print each model in this chapter:

Filterblock	Time_Hr	Time_Min	PLA_Length(m)
FB_Filter_block_simple	2	34	4.58
FB_Filter_collar_compression_tool	0	6	0.26
FB_Filter_collar	0	1	0.03
FB_Filter_F17_slider	0	21	0.58
FB_Filter_F_slider	0	12	0.47
FB_Filter_slider	0	19	0.56
FB_Filter_slot_bottom	0	19	0.38
FB_Filter_slot_stopper	0	8	0.35
FB_Filter_slot_top	0	14	0.33
FB_Infinity_adapter	0	41	0.95
FB_Side_port_separate	0	35	0.99
FB_Side_port_separate_EpiStop	0	51	1.25
FB_Splitter_case	1	32	2.58
FB_Stopper	0	24	0.73
FB_Top_connector	0	42	1.13


## FB\_Filter\_block\_simple




Use tree supports: Branch angle 40. touching baseplate only. Overhang 67 degrees. See picture. The tree gives a solid concentric trunk at the base (better for adhesion). Normal support will result in a hollow square base - less contact area so risk of coming off).




## FB\_Filter\_collar\_compression\_tool




## FB\_Filter\_collar




## FB\_Filter\_F17\_slider




## FB\_Filter\_F\_slider




## FB\_Filter\_slider




## FB\_Filter\_Slot\_bottom




## FB\_Filter\_slot\_stopper




## FB\_Filter\_slot\_top




## FB\_Infinity\_adapter



## FB\_Side\_port\_separate




## FB\_Side\_port\_separate\_EpiStop




Print with supports on: General (not tree) support, overhang 67 degrees, supports enabled 'everywhere'.

## FB\_Splitter\_case




Ensure all supports are off.

## FB\_Stopper



Ensure all supports are off.

## FB\_Top\_connector



## Focus Gears


These are the parts for the stage focus mechanism. The CAD source models for these files are found in the file Focus\_Gears.FCStd.

## Resources


Cura calculates the following resources are required to print each model in this chapter:

Focus Gears	Time_Hr	Time_Min	PLA_Length(m)
FG_Fine_gear	1	38	2.65
FG.Focus_spacer	0	20	0.7
FG_Intermedius	1	7	1.75
FG_Pulley_coarse	1	18	2.03
FG_Pulley	0	21	0.33
FG_Eccentric_Tensioner_Top	0	24	0.69
FG_Eccentric_Tensioner_Bottom	0	19	0.44
FG_Eccentric_Tensioner_Pulley_c_adhesin	0	17	0.22


## FG\_Fine\_gear




## FG\_Focus\_spacer




## FG\_Intermedius




## FG\_Pulley\_coarse




## FG\_Pulley




## FG\_Eccentric\_Tensioner\_Top



## FG\_Eccentric\_Tensioner\_Bottom



## FG\_Eccentric\_Tensioner\_Pulley\_c\_adhesin



## Legs

These are the parts for the various legs / stands options for the microscope. The CAD source models for these files are found in the file Legs.FCStd.


All these should be printed with the 'Plates' custom Cura profile. Only the 'LG\_Feet\_linker' requires supports and special instructions (see below).

## Resources


Cura calculates the following resources are required to print each model in this chapter:

Legs	Time_Hr	Time_Min	PLA_Length(m)
LG_Back_leg_angled	0	36	1.68
LG_Feet_linker	3	18	8.43
LG_Front_legs	3	21	9.34
LG_Hind_extension	0	46	1.64
LG_Short_leg	0	16	0.54

## LG\_Back\_leg\_angled




## LG\_Feet\_linker




Use supports 'touching baseplate only' with 67 degree overhang angle.

Also, you must rotate it in Cura by 90 degrees to the position shown in the figure and ensure it is lain flat agains the build plate because there is a slight angle to it when first loaded (this may be imperceptible by eye but it can cause print failure if not corrected). Both the rotation by 90 degrees and the laying flat can be achieved in one step by selecting the model, going to the rotations menu and click the 'Lay flat' icon – this will do both things for you automatically (tested and works with Cura 4.8 and Cura 4.10.0).


## LG\_Front\_legs



## LG\_Hind\_extension



## LG\_Short\_leg



## Monocular Head

These are the parts for the monocular viewing head of the microscope. Some of the parts here are also used for the binocular and trinocular head modules. The CAD source models for these files are found in the file Monocular.FCStd.

For tubular structures (including the monocular tubes and projection cone) use the 'Flats' profile but modify it thus: use a 'Zig-Zag' infill pattern - it is quicker due to fewer stop-start retracts on the nozzle (16% of total time spent on retractions compared to 40% for the default cubic infill pattern, for the monocular tube model). Also use the 'concentric' 'Top/bottom' pattern.


For other models just use the default 'Flats' profile.

## Resources


Cura calculates the following resources are required to print each model in this chapter:

Monocular Head	Time_Hr	Time_Min	PLA_Length(m)
MN_EM_Block	1	42	2.7
MN_EM_Block_cover	0	31	1.01
MN_Monocular_tube_c_adhesin	2	33	8.04
MN_Monocular_tube_CM_c_adhesin	2	7	6.56
MN_Ocular_cap_170	1	9	2.83
MN_Ocular_cap	0	42	1.65
MN_Ocular_extension_c_adhesin	0	55	2.19
MN_Ocular_extension_CM_c_adhesin	0	38	1.7
MN_Ocular_Extn_CM_170_c_adhesin	1	2	2.62
MN_Ocular_lock_nut	0	16	0.47
MN_Ocular_protective_cap	0	19	1.14
MN_Ocular_tube_protective_cap	0	26	1.14
MN_Projector_cone	2	6	7.91
MN_Aperture_20mm	0	19	0.93
MN_Aperture_46mm	0	7	0.32

## MN\_EM\_Block




## MN\_EM\_Block\_cover




Note the orientation. For some reason, if you try to rotate the mesh by 45 degrees in FreeCAD and export the mesh, the exported mesh needs to be 'laid flat' in Cura because the 45 degree rotation doesn't quite lay it flat.

I have found that it is best to export the mesh from FreeCAD as is, without the 45 degree rotation in FreeCAD, then do that rotation in Cura to position it as shown.


## MN\_Monocular\_tube\_c\_adhesin




## MN\_Monocular\_tube\_CM\_c\_adhesin




## MN\_Ocular\_cap\_170




## MN\_Ocular\_cap




## MN\_Ocular\_extension\_c\_adhesin




## MN\_Ocular\_extension\_CM\_c\_adhesin




## MN\_Ocular\_Extn\_CM\_170\_c\_adhesin




## MN\_Ocular\_lock\_nut




## MN\_Ocular\_protective\_cap




## MN\_Ocular\_tube\_protective\_cap




## MN\_Projector\_cone



## MN\_Aperture\_20mm



## MN\_Aperture\_46mm



# PUMA Control Console


These are the parts for the PUMA Control Console (PCC) used not only to drive the LED light of the microscope but also allowing control of the Z-stage motor and 25x25 mm ST7789 TFT display (the display may be used for either the AR HUD, the SLM or the optional monitor attached to the PCC – only one TFT can be driven with the standard Arduino Nano-based PCC so you can only use one of those options at a time). The CAD source models for these files are found in the file PUMA\_Control.FCStd.

## Resources


Cura calculates the following resources are required to print each model in this chapter:

PUMA Control Console	Time_Hr	Time_Min	PLA_Length(m)
PC_Battery_cover	0	31	1.44
PC_Expansion_port_cover	0	11	0.37
PC_Front_panel	4	21	7.44
PC_Joystick_PCB_clasp	0	18	0.71
PC_Lamp_insulator	0	1	0.02
PC_Left_panel_Base_Skeleton	12	44	24.27
PC_Top_back_panel	8	28	13.42
PC_ArdU_cover	0	6	0.17
PC_Current_knob	0	21	0.71
PC_Monitor_case	3	38	13.73
PC_Monitor_connector	0	24	1.05
PC_Monitor_lens_retainer	0	17	0.89
PC_Monitor_light_shield	0	28	0.79


## PC\_Battery\_cover



## PC\_Expansion\_port\_cover




## PC\_Front\_panel




Use the 'Stamps' Cura profile. Ensure all supports are disabled.


## PC\_Joystick\_PCB\_clasp



## PC\_Lamp\_insulator




## PC\_Left\_panel\_Base\_Skeleton




Ensure all supports are disabled.


## PC\_Top\_back\_panel



Use the 'Stamps' Cura profile and add a support blocker for everything except the back protrusion apertures and the top of the arduino window. The support type should be 'Tree' with a support branch angle of 40.0 and overhang of 67.0 and 'Touching baseplate only'. You should get supports as shown in the following pictures.




## PC\_ArdU\_cover




Use the 'Flats' Cura profile but change the infill to 100%. No supports.

## PC\_Current\_knob




Use the 'Flats' Cura profile but change the infill to 100% and the 'Top/Bottom' pattern to 'Concentric'. No supports.

## PC\_Monitor\_case




Use the 'Flats' Cura profile but change the infill to 10% and 'concentric' pattern of infill. No supports.

## PC\_Monitor\_connector




Use the 'Flats' Cura profile. No supports.

## PC\_Monitor\_lens\_retainer



Use the 'Flats' Cura profile but change the infill to 10% and 'concentric' pattern of infill. Also change and the 'Top/Bottom' pattern to 'Concentric'. No supports.

## PC\_Monitor\_light\_shield



Use the 'Flats' Cura profile but change the infill to 10% and 'concentric' pattern of infill. Also change and the 'Top/Bottom' pattern to 'Concentric'. No supports.

## PUMA Lite


These are the parts for the simple control box use to drive the LED light source of the microscope (but provides no other functions). The CAD source models for these files are found in the file PUMALite.FCStd.

## Resources


Cura calculates the following resources are required to print each model in this chapter:

PUMA_Lite	Time_Hr	Time_Min	PLA_Length(m)
PL_Battery_cover	0	17	0.72
PL_Knob	1	22	1.77
PL_Lamp_insulator	0	1	0.02
PL_Left_Base_Back	8	35	13.52
PL_Top_Front_Right	3	55	7.11

## PL\_Battery\_cover




## PL\_Knob



## PL\_Lamp\_insulator


This is identical to the PC\_Lamp\_insulator of the PUMA Control Console (which see).

## PL\_Left\_Base\_Back




Ensure supports are disabled. Print with the PUMA custom Cura profile 'Stamps'

## PL\_Top\_Front\_Right



Use tree supports with support blocker to prevent the screw hole being supported - as shown in the figure. Support 'everywhere' with overhang 67 degree and tree support branch angle 50 degrees.



# Quick Release Objective Holder


These are the parts for the objective holder that fits onto the focus platform of the microscope stage via a quick-release mechanism. The CAD source models for these files are found in the file QuickRelease\_v2.0.FCStd.

## Resources


Cura calculates the following resources are required to print each model in this chapter:

Quick Release Objective Holder	Time_Hr	Time_Min	PLA_Length(m)
QR_Base_thread	0	43	1.23
QR_C-RMS_Thread	0	12	0.41
QR_Male_C_extn_1mm	1	30	2.47
QR_Male_C_extn_4mm	1	42	2.82
QR_Trainer	1	8	2.5
QR_Trainer_Male	1	25	2.96

## QR\_Base\_thread



## QR\_C-RMS\_Thread




This thread adapter was designed to conform to measurements of similar adapters that are commercially available and made of metal. Thus you can use the metal adapter if you want all-metal threads (which are easier to align with your objective and harder to cross-thread and destroy and wear out).

If using this 3D printed plastic version, the female RMS thread will require training with a metal RMS male thread (such as on an objective lens). Great care is required not to cross the 3D printed plastic thread and damage it at all times but especially during training.

The male C thread is best trained on a metal C thread such as may be found on a C/CS thread extender ring.


## QR\_Male\_C\_extn\_1mm



The female C thread is best trained on a metal C thread such as may be found on a C-mount lens or C/CS thread extender adapter ring. Great care is required to avoid crossing the 3D printed plastic thread and thereby damaging or destroying it.

The outer quick-release male tri-helix thread should be trained using the QR\_Trainer 3D printed tool (see below).


## QR\_Male\_C\_extn\_4mm




The female C thread is best trained on a metal C thread such as may be found on a C-mount lens or C/CS thread extender adapter ring. Great care is required to avoid crossing the 3D printed plastic thread and thereby damaging or destroying it.

The outer quick-release male tri-helix thread should be trained using the QR\_Trainer 3D printed tool (see below).

## QR\_Trainer



## QR\_Trainer\_Male



## Stabiliser

These are the parts for the optional stabiliser exoskeleton bracket. The CAD source models for these files are found in the file Stabiliser.FCStd.


Use the 'Strong01' Cura profile to print the parts.

## Resources

Cura calculates the following resources are required to print each model in this chapter:


Stabiliser	Time_Hr	Time_Min	PLA_Length(m)
SB_S1_Axial	7	29	20.71
SB_S1_Insert_c_adhesin	0	32	0.96
SB_S1_Lateral_strut_Left	2	5	6.26
SB_S1_Lateral_strut_Right	2	6	6.26
SB_S1_Peg	0	16	0.44

## **SB\_S1\_Axial**




Use 'Strong01' profile with concentric top-bottom pattern and with supports 'touching baseplate only'.

## SB\_S1\_Insert\_c\_adhesin




Use 'Strong01' profile with 'lines' top-bottom pattern and with supports 'everywhere' and overhang of 40 degrees.

## SB\_S1\_Lateral\_strut\_Left




The left and right strut models may be printed together (as they are one model in the FreeCAD file) but I recommend they be printed separately as shown here in order to avoid stringing and unnecessary back-forth movements of the print head.

To achieve this, create a cube shape in FreeCAD ('Parts' workbench) and make the cut larger than either one of the struts. Then place the cube over the strut you don't want and perform a subtraction between the struts model and the cube using the shape 'Cut' operation (see the figure below).




Then generate the mesh of the one remaining strut. Once the mesh is made, undo the shape cut, move the cube over the opposite strut and repeat the process to get the other mesh.

## SB\_S1\_Lateral\_strut\_Right



The left and right strut models may be printed together (as they are one model in the FreeCAD file) but I recommend they be printed separately as shown here in order to avoid stringing and unnecessary back-forth movements of the print head. The notes on the previous page show how to do this.

## SB\_S1\_Peg



It may help to set top/bottom layer pattern to 'concentric'

## Stage


These are the parts for the main stage base plate and focus platform of the microscope. The CAD source models for these files are found in the file Stage.FCStd.

## Resources


Cura calculates the following resources are required to print each model in this chapter:

Stage	Time_Hr	Time_Min	PLA_Length(m)
ST_Articulation	0	14	0.53
ST_BasePlate	8	23	25.63
ST_FocusPlate	4	8	11.22

## ST\_Articulation




## ST\_BasePlate



Print with the 'Plates' custom Cura profile.

## ST\_FocusPlate



Print with the 'Plates' custom Cura profile.

## Trinocular Camera Port


This is the part for the trinocular extension that attaches to the advanced filter block and allows a third eyepiece and / or a camera to be attached (in addition to any cameras or eyepieces that may be attached to the main viewing oculars in the ocular head). The CAD source models for these files are found in the file Trinocular\_CP\_v1.FCStd.

## Resources


Cura calculates the following resources are required to print each model in this chapter:

Trinocular_Camera_Port	Time_Hr	Time_Min	PLA_Length(m)
TN_Trinoc_CP_tube	1	37	4.48
TN_Trinoc_CP_CM	1	43	4.47
TN_Ocular_spacer_ring_1p44	0	7	0.28


## TN\_Trinoc\_CP\_tube



## TN\_Trinoc\_CP\_CM



## TN\_Ocular\_spacer\_ring\_1p44



## XY Stabiliser


This is the part for an optional mechanism to make the standard XY mechanical stage calipers more stable in their motion. The CAD source models for this STL file is found in the file XY\_Stabiliser.FCStd.

## Resources

Cura calculates the following resources are required to print each model in this chapter:

XY_Stabiliser	Time_Hr	Time_Min	PLA_Length(m)
XY_Stabiliser	1	26	2.59

## XY\_Stabiliser



Print with the 'XYStabiliser' custom Cura profile.

## Z-Motor


These are the parts for adding a Z-motor to the stage. The CAD source models for these files are found in the file Z\_Motor.FCStd.

## Resources

Cura calculates the following resources are required to print each model in this chapter:


Z-Motor	Time_Hr	Time_Min	PLA_Length(m)
ZM_Motor_attachment	1	1	2.36
ZM_Motor_gear	1	34	2.36
ZM_Z_Limit_Sw_Mount	0	17	0.44
ZM_Z_Probe	1	4	2.04
ZM_Z_Yoke	0	22	0.66

## ZM\_Motor\_attachment




Ensure supports are disabled. It is not possible to properly rotate this in Cura using the method described below (for ZM\_Z\_Limit\_Sw\_Mount) so the mesh is rotated 235 degrees in FreeCAD prior to saving to get the flat parts of the front surface aligned to the X-axis.

## ZM\_Motor\_gear



## ZM\_Z\_Limit\_Sw\_Mount




Ensure supports are disabled. It is difficult to get the flat outer surface exactly in line with the X-axis of the printer. To do this follow these steps:

1. Load into Cura.
2. Rotate it so the flat surface with the screw holes faces the buildplate then use the 'Lay flat' option of Cura to calculate the exact rotation to make that surface flat against the build plate.
3. When done, simply use snap rotation to 90 degrees about the long axis of the model so that the flat surface that was laying on the build plate now faces to the front of the printer.


The resulting position is as shown in the figure.

## ZM\_Z\_Probe



Ensure supports are disabled. It is difficult to get the flat outer surface exactly in line with the X-axis of the printer so a method similar to that described above for 'ZM\_Z\_Limit\_Sw\_Mount' is used. The result is as shown in the picture.

## ZM\_Z\_Yoke



# CNC-Stage

These are the parts for making the CNC precision motorised XYZ stage. The CAD source models for these files are found in the file CNC\_Stage.FCStd.

All models are printed with the basic Cura profile ‘Strong02’ (10 wall layers and 10 layers top and bottom) and most are printed with an infill of 30%. Details of any modifications to this and any supports / support blockers will be provided with each model.

## Location of the models

### File: CNC\_Stage.FCStd

#### *Group: PSU*

- PSU\_CableFixture\_front
- PSU\_DC\_CableCover
- PSU\_Mains\_CableCover
- PSU\_Peg\_Leg

#### *Group: Arduino*

- Ardu\_bus\_bar\_mount
- Ardu\_Bus\_Clip\_set
- Ardu\_M3\_Washer

#### *Group: X\_Axis*

- X\_Lim\_Sarcoph\_L
- X\_Lim\_Sarcoph\_R
- X\_NEMA11\_Bracket
- X\_NEMA11\_Counter\_weight
- X\_NEMA11\_CW\_Hopper
- X\_Probe
- X\_Probe\_Housing\_L
- X\_Probe\_Housing\_R
- X\_Probe\_shim
- XY\_NEMA11\_Undertable\_Support

#### *Group: Y\_Axis*

- Y\_Lim\_Sarcoph\_Susp\_B
- Y\_Lim\_Sarcoph\_Susp\_F
- Y\_NEMA11\_Bracket

Y\_Probe\_Susp\_Nemesh  
Y\_Suspension\_platform\_B  
Y\_Suspension\_platform\_F

*Group: Z\_Axis*

Z\_Coupler  
Z\_Lim\_Attachment  
Z\_Lim\_Back\_baffle  
Z\_Lim\_Probe  
Z\_M4\_small\_washer  
Z\_Mount  
Z\_N11\_Syringe\_body  
Z\_N11\_Syringe\_cover  
Z\_N11\_Syringe\_Gearbox\_plunger  
Z\_N11\_Syringe\_mount

*Group: Sample\_holder*

Sam\_Counter\_weight  
Sam\_Pincer\_Extn1\_L  
Sam\_Pincer\_Extn1\_R  
Sam\_Pincer\_Extn2\_L  
Sam\_Pincer\_Extn2\_R  
Sam\_Presser\_wings  
Sam\_Spring\_arm\_assy  
Sam\_Stage\_Top

*Group: BaseBoard*

BB\_Jig\_Back\_Left  
BB\_Jig\_Back\_Measure  
BB\_Jig\_Back\_Right  
BB\_Jig\_Front\_Left  
BB\_Jig\_Front\_Measure  
BB\_Jig\_Front\_Right  
BB\_Jig\_Left\_Measure  
BB\_Jig\_Right\_Measure  
BB\_Jig\_Undercarriage  
BB\_Jig\_XYTable  
BB\_Jig\_Y\_Lim

*Group: Frame*

- F\_Ardu\_strut\_BL
- F\_Ardu\_strut\_BR
- F\_Ardu\_strut\_TL
- F\_Ardu\_strut\_TR
- F\_DM3\_Mount
- F\_OM\_QR\_Plate
- F\_Overmount
- F UM Opt\_gripper
- F UM Opt\_Gripper\_washer
- F UM Opt\_to\_QR\_F
- F\_Undercarriage\_platform

*Group: Shields*

- SD\_2020\_End\_cap
- SD\_2020\_Protector
- SD\_Aft\_clip
- SD\_Aft\_Hole\_cover
- SD\_Aft\_Left
- SD\_Aft\_Mid
- SD\_Aft\_Right
- SD\_Aft\_Rivet\_block
- SD\_DIN\_End\_Protector
- SD\_DIN\_Foot\_Protector
- SD\_Dorsal\_Key\_head
- SD\_Dorsal\_Key\_key
- SD\_Dorsal\_Key\_main
- SD\_Dorsal\_Key\_wedge
- SD\_Dorsal\_L
- SD\_Dorsal\_M
- SD\_Dorsal\_R
- SD\_Dorsal\_Rivet\_plate
- SD\_Dorsal\_suture
- SD\_Forward\_L
- SD\_Forward\_R
- SD\_Forward\_Rivet\_handle\_1

SD\_Forward\_Rivet\_handle\_2  
 SD\_Forward\_Rivet\_lock  
 SD\_Forward\_Rivet\_tool  
 SD\_Forward\_suture  
 SD\_Lateral\_Left  
 SD\_Lateral\_Right  
 SD\_Letters\_CNC  
 SD\_Letters\_PUMA  
 SD\_T\_nut\_positioner  
 SD\_XMotor\_cover\_L  
 SD\_XMotor\_cover\_R

## Resources

Cura 4.13.1 calculates the following resources are required to print each model in this chapter:


CNC-Stage	Time_Hr	Time_Min	PLA_Length(m)
CN_PSU_CableFixture_front	9	1	29.11
CN_PSU_DC_CableCover	0	32	1.66
CN_PSU_Mains_CableCover	1	1	2.99
CN_PSU_Peg_Leg	0	9	0.3
CN_Ardubus_bar_mount	7	29	20.74
CN_Ardubus_Clip_set	4	13	6.1
CN_Ardum3_Washer	0	1	0.02
CN_X_Lim_Sarcoph_L	2	35	7.11
CN_X_Lim_Sarcoph_R	2	30	7.12
CN_X_NEMA11_Bracket	9	41	21.55
CN_X_NEMA11_Counter_weight	10	32	25.04
CN_X_NEMA11_CW_Hopper	3	27	11.69
CN_X_Probe	0	26	1.04
CN_X_Probe_Housing_L	0	53	2.77
CN_X_Probe_Housing_R	0	54	2.77
CN_X_Probe_shim	0	2	0.1

CN_XY_NEMA11_Undertable_Support	0	34	1.72
CN_Y_Lim_Sarcoph_Susp_B	4	36	12.35
CN_Y_Lim_Sarcoph_Susp_F	4	36	12.35
CN_Y_NEMA11_Bracket	9	7	20.82
CN_Y_Probe_Susp_Nemesh	3	53	9.31
CN_Y_Suspension_platform_B	0	36	1.97
CN_Y_Suspension_platform_F	0	36	1.97
CN_Z_Coupler	0	37	2.2
CN_Z_Lim_Attachment	2	23	6.44
CN_Z_Lim_Back_baffle	0	33	0.62
CN_Z_Lim_Probe	0	16	0.6
CN_Z_M4_small_washer	0	0.4	0.02
CN_Z_Mount	10	17	28.73
CN_Z_N11_Syringe_body	6	11	14.84
CN_Z_N11_Syringe_cover	0	42	1.66
CN_Z_N11_Syringe_Gearbox_plunger	2	45	7.09
CN_Z_N11_Syringe_mount	8	26	23.77
CN_Sam_Counter_weight	1	6	3.15
CN_Sam_Pincer_Extn1_L	0	52	2.3
CN_Sam_Pincer_Extn1_R	0	52	2.3
CN_Sam_Pincer_Extn2_L	1	24	3.92
CN_Sam_Pincer_Extn2_R	1	24	3.94
CN_Sam_Presser_wings	0	53	2.47
CN_Sam_Spring_arm_assy	6	1	15.43
CN_Sam_Stage_Top	11	0	28.98
CN_BB_Jig_Back_Left	5	10	12.8
CN_BB_Jig_Back_Measure	1	13	3.02
CN_BB_Jig_Back_Right	5	10	12.8
CN_BB_Jig_Front_Left	4	57	12.47
CN_BB_Jig_Front_Measure	0	17	0.87
CN_BB_Jig_Front_Right	4	57	12.47
CN_BB_Jig_Left_Measure	1	22	3.6
CN_BB_Jig_Right_Measure	1	22	3.6


CN_BB_Jig_Undercarriage	4	49	12.12
CN_BB_Jig_XYTable	5	33	13.98
CN_BB_Jig_Y_Lim	0	38	2.09
CN_F_Ardु_strut_BL	1	49	5.8
CN_F_Ardु_strut_BR	1	49	5.8
CN_F_Ardु_strut_TL	1	49	5.8
CN_F_Ardु_strut_TR	1	49	5.8
CN_F_DM3_Mount	2	12	6.32
CN_F_OM_QR_Plate	9	43	27.37
CN_F_Overmount	9	33	26.48
CN_F_UM_Opt_gripper	1	42	4.63
CN_F_UM_Opt_Gripper_washer	0	3	0.16
CN_F_UM_Opt_to_QR_F	2	40	7.81
CN_F_Undercarriage_platform	8	36	23.99
CN_SD_2020_End_cap	0	13	0.46
CN_SD_2020_Protector	0	23	0.9
CN_SD_Aft_clip	0	5	0.07
CN_SD_Aft_Hole_cover	0	34	1.84
CN_SD_Aft_Left	16	13	67.98
CN_SD_Aft_Mid	15	2	54.76
CN_SD_Aft_Right	16	12	67.98
CN_SD_Aft_Rivet_block	0	9	0.39
CN_SD_DIN_End_Protector	0	20	0.86
CN_SD_DIN_Foot_Protector	0	48	2.3
CN_SD_Dorsal_Key_head	0	10	0.38
CN_SD_Dorsal_Key_key	0	7	0.35
CN_SD_Dorsal_Key_main	0	4	0.1
CN_SD_Dorsal_Key_wedge	0	0.4	0.01
CN_SD_Dorsal_L	15	51	47.81
CN_SD_Dorsal_M	8	40	29.69
CN_SD_Dorsal_R	15	56	48.14
CN_SD_Dorsal_Rivet_plate	0	40	2.1
CN_SD_Dorsal_suture	0	16	0.68

CN_SD_Forward_L	3	50	11.71
CN_SD_Forward_R	3	52	11.72
CN_SD_Forward_Rivet_handle_1	0	4	0.1
CN_SD_Forward_Rivet_handle_2	0	4	0.1
CN_SD_Forward_Rivet_lock	0	9	0.28
CN_SD_Forward_Rivet_tool	0	3	0.12
CN_SD_Forward_suture	0	40	1.86
CN_SD_Lateral_Left	22	47	78.26
CN_SD_Lateral_Right	22	30	78.23
CN_SD_Letters_CNC	0	11	0.52
CN_SD_Letters_PUMA	0	15	0.74
CN_SD_T_nut_positioner	0	4	0.07
CN_SD_XMotor_cover_L	9	2	36.4
CN_SD_XMotor_cover_R	9	2	36.4


## CN\_PSU\_CableFixture\_front




## CN\_PSU\_DC\_CableCover




## CN\_PSU\_Mains\_CableCover




## CN\_PSU\_Peg\_Leg



## CN\_Ardु\_bus\_bar\_mount




## CN\_Ardu\_Bus\_Clip\_set




Note: Printing these as a set is convenient, time-saving and allows for better print-bed adhesion but you could print the individual clips in this set separately if you want to. They must be orientated on the print bed on their 'side' as shown in this image, even if only printing one alone.

## CN\_Ardu\_M3\_Washer




Use 'Concentric' for Top/bottom pattern.

## CN\_X\_Lim\_Sarcoph\_L




This uses tree supports 'everywhere' as shown in the top inset.

## CN\_X\_Lim\_Sarcoph\_R




This also uses tree supports (as for CN\_X\_Lim\_Sarcoph\_L) but, in addition, there is a support blocker stopping the supports going in the M2 holes and limit switch alcove, as shown.

## CN\_X\_NEMA11\_Bracket




This must be printed on its 'side', as shown, in order to achieve the correct strength in the correct directions. Note the position of the support blockers. This uses tree supports 'everywhere' with a support overhang angle of 60 degrees. This also uses the concentric pattern for top and bottom layers and zig-zag infill pattern with an infill density of 25%.

## CN\_X\_NEMA11\_Counter\_weight




This must be printed on its 'side', as shown, in order to achieve the correct strength in the correct directions. Note the position of the support blockers. This uses normal supports as shown.

## CN\_X\_NEMA11\_CW\_Hopper




This uses normal supports 'touching buildplate only' with a 67 degree overhang limit.

## CN\_X\_Probe




This uses normal supports upto 67 degrees 'touching buildplate only'. Note the support blocker in the holes.

## CN\_X\_Probe\_Housing\_L




This uses normal supports upto 67 degrees 'touching buildplate only'.

## CN\_X\_Probe\_Housing\_R




This uses normal supports upto 67 degrees 'touching buildplate only'.

## CN\_X\_Probe\_shim




## CN\_XY\_NEMA11\_Undertable\_Support



Use a 'Top/bottom Pattern' of 'Concentric'.


## CN\_Y\_Lim\_Sarcoph\_Susp\_B



This is best printed together with the CN\_Y\_Lim\_Sarcoph\_Susp\_F (which see for illustration of print bed orientation). Use tree supports and support blockers in the M2 holes just as described for the CN\_X\_Lim\_Sarcoph\_R model (which see for details).

The figure illustrates the orientation of the models for fitting on the stage base board. See the section on the CN\_Y\_Lim\_Sarcoph\_Susp\_F for details.


## CN\_Y\_Lim\_Sarcoph\_Susp\_F



This is best printed together with the CN\_Y\_Lim\_Sarcoph\_Susp\_B as shown. Use tree supports and support blockers in the M2 holes just as described for the CN\_X\_Lim\_Sarcoph\_R model (which see for details).


Note on orientation of the printed models: The suspension side-car has a rounded corner adjacent to its sarcophagus at the back and a squared off corner connection to its sarcophagus at the front. This is illustrated in the CN\_Y\_Lim\_Sarcoph\_Susp\_B section (which see). The front is where the XY table is situated.

## CN\_Y\_NEMA11\_Bracket



This is printed on its 'side' with support blockers as shown and using tree supports. This is similar to the CN\_X\_NEMA11\_Bracket (which see for more detail).

## CN\_Y\_Probe\_Susp\_Nemesh




This is printed with normal supports ‘everywhere’ with a 67 degree overhand angle and with support blockers as shown.

## CN\_Y\_Suspension\_platform\_B

This is the same as the CN\_Y\_Suspension\_platform\_F model and the two are best printed together as shown in the illustration below.


## CN\_Y\_Suspension\_platform\_F



This is best printed with the CN\_Y\_Suspension\_platform\_B, as shown. The two models are actually identical, separation into 'F' and 'B' types is purely for admin purposes.


The orientation may not seem optimal (and it is not optimal for part strength in the thin end-parts) but it is optimal for producing a flat surface for the suspension rollers to traverse (the surface that is applied to the print bed) and hence for producing a 'smooth ride' for the precision mechanism. The thin end parts will not be under strain when fully fitted into their suspension side-cars so optimal strength is not needed here and so strength can be sacrificed for the benefit of precision – hence the preferred build plate orientation is as shown.

## CN\_Z\_Coupler



Use 'concentric' for both the Top/bottom pattern and the infill pattern. No supports are used.


## CN\_Z\_Lim\_Attachment



Use tree supports and support blockers as shown.


Note that I am also printing the CN\_Z\_Lim\_Back\_baffle model in here to save time. You can print it separately if you wish (see CN\_Z\_Lim\_Back\_baffle).

## CN\_Z\_Lim\_Back\_baffle




Some form of supports must be used. Note that you could either print it separately as shown here or you could print it altogether with the CN\_Z\_Lim\_Attachment model – as is illustrated in the figure in the section for that model.

## CN\_Z\_Lim\_Probe




## CN\_Z\_M4\_small\_washer



Use 'concentric' for 'Top/bottom pattern'.


Might as well print multiple of these all at once because it is very quick to print and you will need at least 4.

## CN\_Z\_Mount




Use 'Normal' supports 'touching buildplate only' with an overhang angle of 67 degrees.


## CN\_Z\_N11\_Syringe\_body




## CN\_Z\_N11\_Syringe\_cover



## CN\_Z\_N11\_Syringe\_Gearbox\_plunger




## CN\_Z\_N11\_Syringe\_mount



Use 'Normal' supports 'touching buildplate only' with an overhang angle of 67 degrees.

## CN\_Sam\_Counter\_weight




Use 'Normal' supports 'everywhere' with a 67 degree overhang angle.

## **CN\_Sam\_Pincer\_Extn1\_L**

This will usually be printed together with the R version, as shown below.

## **CN\_Sam\_Pincer\_Extn1\_R**

This will usually be printed with the L version, as shown below.




## **CN\_Sam\_Pincer\_Extn2\_L**


This will usually be printed together with the R version, as shown below.

## **CN\_Sam\_Pincer\_Extn2\_R**

This will usually be printed together with the L version, as shown below.




## CN\_Sam\_Presser\_wings




Use 'Normal' supports at 67 degrees and 'Touching buildplate only'


## CN\_Sam\_Spring\_arm\_assy



## CN\_Sam\_Stage\_Top




## CN\_BB\_Jig\_Back\_Left



This shows the arrangement for printing both the left and right 'Jig\_Back' components together.


## CN\_BB\_Jig\_Back\_Measure



## **CN\_BB\_Jig\_Back\_Right**


See the section for 'CN\_BB\_Jig\_Back\_Left'

## CN\_BB\_Jig\_Front\_Left



This shows the arrangement for printing both the left and right 'Jig\_Front' components together.


## CN\_BB\_Jig\_Front\_Measure



## **CN\_BB\_Jig\_Front\_Right**

See the section for 'CN\_BB\_Jig\_Front\_Left'

## **CN\_BB\_Jig\_Left\_Measure**




This shows the arrangement for printing both the left and right components together.


## **CN\_BB\_Jig\_Right\_Measure**

See the section on 'CN\_BB\_Jig\_Left\_Measure'


## CN\_BB\_Jig\_Undercarriage



## CN\_BB\_Jig\_XYTable



## CN\_BB\_Jig\_Y\_Lim



## **CN\_F\_Ardu\_strut\_BL**

These are best printed all together as shown in the section on 'CN\_F\_Ardu\_strut\_TR'.


## **CN\_F\_Ardu\_strut\_BR**

These are best printed all together as shown in the section on 'CN\_F\_Ardu\_strut\_TR'.

## **CN\_F\_Ardu\_strut\_TL**

These are best printed all together as shown in the section on 'CN\_F\_Ardu\_strut\_TR'.

## **CN\_F\_Ardu\_strut\_TR**




Note that this shows all the Ardu struts (not just the TR) because it is recommended to print them all together as shown. Note also that each of these struts has a unique shape, none of them are duplicates or mirror images of any other (so don't try to just load one then do 'multiply-model').

## CN\_F\_DM3\_Mount




Note that this shows 3 copies of the DM3\_Mount printed all together. This is recommended if you are using one DM320T (or DM320T form-factor compatible) digital microstepper driver for each of the three axes. If you are using a different form factor driver or CNC motherboard then these models will not be needed and it will be up to you to design appropriate mounts for your electronics.

## CN\_F\_OM\_QR\_Plate




This uses 'Normal' supports 'touching buildplate only' with an overhang of 67 degrees.


## CN\_F\_Overmount



## CN\_F\_UM\_Opt\_gripper




## CN\_F UM Opt Gripper washer




Use 'Concentric' for 'Top/Bottom Pattern'.

## CN\_F\_UM\_Opt\_to\_QR\_F




Use 'Concentric' for 'Top/Bottom Pattern'.

## CN\_F\_Undercarriage\_platform




This uses 'Normal' supports 'touching buildplate only' with an overhang of 67 degrees.

## CN\_SD\_2020\_End\_cap




## CN\_SD\_2020\_Protector



## **CN\_SD\_Aft\_clip**

It is recommended to print two of these together with two of the 'CN\_SD\_Aft\_Hole\_cover' model (see below).

## **CN\_SD\_Aft\_Hole\_cover**




The image shows two of this model together with two of the " " model all printed together.


## **CN\_SD\_Aft\_Left**

It is recommended to print this together with the 'CN\_SD\_Aft\_Right' model. See the section on 'CN\_SD\_Aft\_Right' for the build-plate orientation image and further information.

## CN\_SD\_Aft\_Mid




## CN\_SD\_Aft\_Right




The image shows both this model and the 'CN\_SD\_Aft\_Left' model on the build plate together (which is the recommended way to print these in order to save time and energy). They are separated by just 1 mm. No supports are needed.

## CN\_SD\_Aft\_Rivet\_block




The illustration shows 14 copies of this model laid out for printing in one job (you will need 14 of them).


## CN\_SD\_DIN\_End\_Protector



## CN\_SD\_DIN\_Foot\_Protector




## CN\_SD\_Dorsal\_Key\_head



This is best printed as a set of 3 models 'CN\_SD\_Dorsal\_Key\_main', 'CN\_SD\_Dorsal\_Key\_wedge' and 'CN\_SD\_Dorsal\_Key\_head' as shown. Note that the 'CN\_SD\_Dorsal\_Key\_main' model is orientated with the flattened part of its shaft facing the print bed.

## CN\_SD\_Dorsal\_Key\_key




## **CN\_SD\_Dorsal\_Key\_main**

This is best printed as a set of 3 models ‘CN\_SD\_Dorsal\_Key\_main’, ‘CN\_SD\_Dorsal\_Key\_wedge’ and ‘CN\_SD\_Dorsal\_Key\_head’ as shown in the illustration in the section on ‘CN\_SD\_Dorsal\_Key\_head’ (which, see).

## **CN\_SD\_Dorsal\_Key\_wedge**


This is best printed as a set of 3 models ‘CN\_SD\_Dorsal\_Key\_main’, ‘CN\_SD\_Dorsal\_Key\_wedge’ and ‘CN\_SD\_Dorsal\_Key\_head’ as shown in the illustration in the section on ‘CN\_SD\_Dorsal\_Key\_head’ (which, see).

## CN\_SD\_Dorsal\_L




Use 'Tree' supports 'touching buildplate only'.

## CN\_SD\_Dorsal\_M




## CN\_SD\_Dorsal\_R




Use 'Tree' supports 'touching buildplate only'.


## CN\_SD\_Dorsal\_Rivet\_plate



## CN\_SD\_Dorsal\_suture



## CN\_SD\_Forward\_L




The left and right forward shields may be printed together with the forward suture, as shown. You will need to enable 'Normal' supports at an overhang angle of 67 degrees for the shields (the suture does not need this).

## CN\_SD\_Forward\_R

See the section on 'CN\_SD\_Forward\_L'.

## CN\_SD\_Forward\_Rivet\_handle\_1



The 'rivet-lock' system for the forward shields are best printed as a set, as shown in the illustration. The system consists in one 'Tool' and 10 lots of handled rivet locks, each rivet lock paired with one handle but 5 of the handles are of type 1 (with the flattened part of the aperture closest to the handle) and the other 5 handles are of type 2 (the flat part of the aperture faces away from the handle).

## CN\_SD\_Forward\_Rivet\_handle\_2

See 'CN\_SD\_Forward\_Rivet\_handle\_1'.

## CN\_SD\_Forward\_Rivet\_lock

See 'CN\_SD\_Forward\_Rivet\_handle\_1'.


## CN\_SD\_Forward\_Rivet\_tool

See 'CN\_SD\_Forward\_Rivet\_handle\_1'.

## **CN\_SD\_Forward\_suture**


See the section on 'CN\_SD\_Forward\_L'.

## CN\_SD\_Lateral\_Left




Use 'Normal' supports 'touching buildplate only' with a 67 degree overhang.

## CN\_SD\_Lateral\_Right




Use 'Normal' supports 'touching buildplate only' with a 67 degree overhang.

## CN\_SD\_Letters\_CNC




Submerge the model 2.8 mm below the build plate surface before slicing (as shown) and use the 'Concentric' value for 'Top/Bottom pattern'.

## CN\_SD\_Letters\_PUMA




Submerge the model 2.8 mm below the build plate surface before slicing (as shown) and use the 'Concentric' value for 'Top/Bottom pattern'.

## CN\_SD\_T\_nut\_positioner



These are very small and quick to print and you will need 6 of them so it makes sense to print them all together, as shown.

## CN\_SD\_XMotor\_cover\_L



The figure shows both the 'CN\_SD\_Xmotor\_cover\_L' and 'CN\_SD\_Xmotor\_cover\_R' models positioned on the build plate to be printed together for convenience.

Note that the two models are not simply the same model rotated but one is a mirror image of the other and some features (notably two of the articulating keys at the top) will be in the wrong position if you try to duplicate one of the models instead of printing left and right models – so don't do that.

## CN\_SD\_Xmotor\_cover\_R

See the section on 'CN\_SD\_Xmotor\_cover\_L', above.